

# Notice Inviting e-Tender

# West Bengal Medical Services Corporation Limited Swasthya Sathi GN-29, Salt Lake, Sector-V

Kolkata-700091

Phone No (033) 40340307/320 E mail: procurement@wbmsc.gov.in

Supply and Commissioning of Medical equipment for the District Integrated Public Health Laboratory

(Discipline-Microbiology)

(Submission of Bid through online)

Bid Reference No.: WBMSCL/NIT-156/2023

Dated-04.04.2023

# AMENDMENT-II Revised Technical Specification

<u>ltem No-1</u>

# **Biosafety cabinet Class II A2**

- 1. Vertical laminar flow air to create ISO-Class III work area
- 2. Single piece wall of SS, preferably SS- 304
- 3. Single piece work tray of SS, preferably SS-304.
- 4. Side walls are surrounded by negative pressure.
- 5. LUPA filter (efficacy 99.999%. Arrests particles \$ 0.3micron) 2 (both down flow and exhaust).
- 6. UV light arrangement with a UV meter.
- 7. Pressure gauge (dial type).
- 8. Size :  $WxDx H = \frac{4 \text{ ft. } x 2 \text{ ft. } x 2 \text{ ft. } (approx)}{4 \text{ ft. } x 2 \text{ ft. } (approx)}$
- 9. Air flow velocity: Inflow 0.5 M/S, Outflow 0.3 M/S 1 (). Heavy Duty Castor wheels
- 10.Particle counter, velocity, calibration certificate (in-built) should be provided.
- 11.Certification: CE / EC / UK Cert / US FDA / BIS

# <u>ltem No-2</u>

# Centrifuge table top (swing out) with 8 tubes upto 6000 RPM

1) Rotor capacity 12 (or more) x 5 ml swing out.

- 2) Maximum speed should be 6000rpm
- 3) Power supply 230 volt +/-10v, 50 Hz AC
- 4) Should have automatic rotor identification
- 5) Should have stainless steel centrifuae chamber

6) Brushless induction motor with variable frequency drive

7) Microprocessor controller with digital display

8) Stable speed output even under unstable voltage condition

9) Smooth and soft start, low temperature rise

### 10) LED display of speed

11) Digital countdown timer and continuous run

- 12) Safety lid interlock to prevent lids opening during centrifugation
- 13) Imbalance detection and centrifugation stop with display of error
- 14) Dynamic break for quick declaration
- 15) Motor over head protection
- 16) Gas hinge to prevent door falling
- 17) Emergency lid lock release
- 18) Should have two (02) year warranty
- 19) Certification: CE / EC / UK Cert / US FDA / BIS

## <u>Item No-3</u>

## Electronic Balance up to 3 (three) decimal places

- 1. Type of Electronic Weighing Machine: electronic weighing scale
- 2. Range: 1 mg to 500 mg
- 3. Readability: 0.1 mg
- 4. Display: LED/ LCD
- 5. Number of Display: single
- 6. Speed of measurement: 1.0 or higher
- 7. Internal calibration: Yes/No
- 8. Tare provision: Manual / Automatic
- 9. Connectivity to Computer: Yes/ No

10.Number of USB ports: 2, NA

- 11.Power Requirement: Mains operated
- 12.Power Supply: Single Phase
- 13.Type of Sensing Element: Load Cell/EMFC
- 14.Warranty in Years: 2
- 15.Pan/ Platform Shape: rectangular/ Square
- 16.Pan/ Platform Size: 290 mm x 240 mm, 300 mm x 300 mm
- 17.Warranty: 2 years
- 18.Legal meteorological license
- 19.Certification: CE / EC / UK Cert / US FDA / BIS

# Item No-4 Hot plate for culture media preparation

- 1. Stainless Steel/ MS power coding Body
- 2. Shape: Rectangular / circular
- 3. Glass ceramic area : 235 mm or more length/ Dia.
- 4. Digital PID controller
- 5. Temperature control: 200 degree Celsius
- 6. Surface Temperature: 250 degree Celsius
- 7. Power requirement: 220V single phase, power consumption ≤1 KW
- 8. Chemical resistance & temperature uniformity
- 9. Certification: CE / EC / UK Cert / US FDA / BIS

#### <u>ltem No-5</u>

#### Microcentrifuge Machine ( up to 16,000 RPM)

- 1. It should be a digital microcentrifuge designed for small sample volume.
- 2. It should have sturdy and durable metal housing (SS 304 grade).
- 3. It should have auto-balancing system.
- 4. Rotor speed should be up to 16000 RPM.
- 5. Should accommodate at least 12 place x 2.0/0.5mL rotor with screw-on lid microtubes
- 6. It should have safety interlock that can prevent the lid from being opened during Operation.
- 7. Overheating of motor should be prevented.
- 8. It should have a precision digital timer that can be set from 0 to 15 minutes at 1 min. increment.
- 9. It should have LED digital display that enables user to monitor the present speed or the remaining time.
- 10. It should have automatic brake system for rapid deceleration.
- 11. The lid can be opened automatically/ through switch once the rotor stops spinning.
- 12. It should accommodate micro or microhematocrit rotors.
- 13. Should be supplied with lid cover (01)
- 14. Certification: CE / EC / UK Cert / US FDA / BIS.
- 15. It should have warranty for 2 years.

#### <u>ltem No-6</u>

#### PCR Workstation

- 1. Work area size 4 ft W x 2ft D x 2ft H.
- 2. Vertical laminar flow design.
- 3. Single piece wall of SS, preferably SS- 304
- 4. Single piece work tray of SS, preferably SS-304.
- 5. Average down flow 60fpm +/- 5fpm at 10cm above the work access opening.

- 6. Down flow velocity should be uniform across the work area and the variance should not be more than 20% from the mean.
- 7. Continuous monitoring and digital display of down flow velocity desirable.
- 8. 254 nm UV lights, one on the work area and two behind the HEPA filter.
- 9. Safety interlock for work area UV light.
- 10. Should have hour counters for HEPA filter and UV lights.
- 11. Mini pleat HEPA filter with a typical efficiency of 99.99% at 0.3 microns or better.
- 12. Work access should be upward sliding type with variable height adjustment.
- 13. Front and sides should be tempered glass.
- 14.2 nos. of 3 pin, 5A electrical sockets with switches should be provided on the back wall.
- 15. Should be provided with support stand with castor wheels and lock.
- 16. Certification: CE / EC / UK Cert / US FDA / BIS

## <u>ltem No-7</u>

### Micro incinerator for inoculating loops & needles

- 1. Max temperature : 800° c or more
- 2. Chamber: ceramic
- 3. Power requirement: 220 V AC Single phase,
- 4. Average cycle time (minutes): 0.1 to 0.2
- 5. Warranty 02 years.
- 6. Certification: CE / EC / UK Cert / US FDA / BIS

### <u>ltem No-8</u>

# Autoclave for sterilization / disinfection – 100 Litre

- 1. Sterilizer chamber capacity in L (usable volume) 100 liter
- 2. Sterilizer Chamber Type: Circular
- 3. Working Pressure of Chamber: 16 pound per sqare inch
- 4. Type of Sterilizer Chamber Door Heavy duty SS 304, hinged door
- 5. Door sealing: By elastomeric rubber gasket suitable to withstand temperature up to 140 degree C & pressure up to 20-30 psi
- 6. Audio- visual Alarm facility available for notifying: Temperature sensor with range 100°C to 300°C, unit should be provided with a low water alarm system as well as steam release valves and safety valves, Cut off automatically when the autoclave is dry
- 7. Working Temperature: 121 degree Celsius
- 8. Range of Pressure Gauge: 30 pound per sqare inch
- 9. Accessories: Water level indicator, Water inlet & outlet valves, Automatic pressure control switch, Mechanical timer, Manual Water Filling & Removal
- 10. The unit must have heater fitted at the bottom and with capacity of more than 3 KW 11. Certification: CE / EC / UK Cert / US FDA / BIS.